

Keep Your Lab Flowing:

Routine Maintenance Supplies for

Waters ACQUITY LCs

Replace Every Month

Bluestem Glass Solvent Filter

- Restek Bluestem glass solvent filter provides clean mobile phase to extend the life of columns and pump seals.
- 15 µm borosilicate glass frit sits lower than conventional glass filters to draw more mobile phase from each bottle.
- Blue filter stem allows instant visual confirmation of upright filter orientation.
- Connects to standard 1/8" OD (3.2 mm) PTFE tubing using your existing frit adaptor (also sold separately as cat.# 26392).

Prevent the particulates and microbial growth in your LC solvents from entering your instrument with the new Restek Bluestem glass solvent filter.

Glass Solvent Filter, 15 µm frit	5041-2168	ea.	26431
Frit Adaptor, PTFE	5062-8517	4-pk.	26392



Solvent Bottle Filter for Waters ACQUITY Systems

Solvent Bottle Filters	ACQUITY	700003616	7-pk.	25954
------------------------	---------	-----------	-------	-------



Replace Every 3 Months

Head Plunger Seal Kits for Waters HPLC Systems

Head Plunger Seal Kit	ACQUITY, nanoACQUITY	700002599	kit	26428
-----------------------	----------------------	-----------	-----	-------



Replace Every 6 Months

Preventive Maintenance (PM) Kits for Waters HPLC & ACQUITY UPLC Systems

ACQUITY BSM Pump Kit

Includes: in-line filter assembly, stainless steel frit; tube assembly, transducer to check valve; tube assembly, SSV to in-line filter; primary check valve (2); wash seal, float flanged, (2); head plunger seal kit (2); sapphire plungers, (2); air filter, pump; air filter, pump handle; mixer assembly, 50 µl; accumulator check valve, double ball & seat (2); solvent bottle filter, stainless steel (7); pump O-ring, (PTFE)

ACQUITY
UPLC BSM

201000173

kit

25797



25797

ACQUITY I2V BSM Pump Kit


Includes: tube assembly, SSV to I2v; transducer to check valve, tube assembly, I2v; wash seal, float flanged (2); head plunger seal kit (2); sapphire plungers (2); air filter, pump; air filter, pump handle; mixer assembly, 50 µl; filter frit cartridge, stainless steel; accumulator check valve, double ball & seat (2); check valve cartridge (2); solvent bottle filter, stainless steel (7); pump O-ring, (PTFE)

ACQUITY
UPLC I2V BSM

201000197

kitv

25798



25798

ACQUITY H-Class QSM Pump Kit

Includes: solvent bottle filter, stainless steel (5); 20 micron frit holder assembly; tube assembly, transducer to check valve; wash seal, float flanged (2); head plunger seal kit (2); sapphire plungers (2); check valve, double ball & seat (2); I2 check valve cartridge; air filter, door; pump O-ring, (PTFE); mixer assembly, 100 µl

ACQUITY
H-CLASS QSM

201000233

kit

25799



25799

Primary Check Valve for Waters ACQUITY and nanoACQUITY Systems

Primary Check Valve	ACQUITY, nanoACQUITY	700002596	2-pk.	25955
---------------------	----------------------	-----------	-------	-------



25955

Accumulator Check Valve (Double Ball & Seat) for Waters ACQUITY and nanoACQUITY Systems

Accumulator Check Valve (Double Ball & Seat)	ACQUITY, nanoACQUITY	700002968	2-pk.	25956
--	----------------------	-----------	-------	-------



25956

Replace Every 12 Months

Sapphire Plunger

for Waters ACQUITY and nanoACQUITY Systems

Sapphire Plunger	ACQUITY, nanoACQUITY	700002600	2-pk.	25958
------------------	----------------------	-----------	-------	-------



25958

Pump O-Ring

for Waters ACQUITY and nanoACQUITY Systems

Pump O-Ring	ACQUITY, nanoACQUITY	WAT076152	ea.	25951
-------------	----------------------	-----------	-----	-------



25951

Needle Assembly

for Waters ACQUITY Systems

Needle Assembly, 0.010 mm ID	ACQUITY	700002644	ea.	25942
------------------------------	---------	-----------	-----	-------



25942

Needle Kit

for Waters ACQUITY Systems

Needle Kit, 15 µL	ACQUITY H-class	700005215	kit	25777
-------------------	-----------------	-----------	-----	-------



25777

Don't Forget!

Lamps for Waters Detectors

Lamp, UV, Long Life (2,000 hours)	ACQUITY PDA/ TUV 2489/2998	201000281	ea.	25775
-----------------------------------	-------------------------------	-----------	-----	-------



25775

This is just the beginning!

There are more RM supplies for your Waters, as well as parts for Agilent, PerkinElmer, Shimadzu, Thermo, and other LCs, at www.restek.com/LCacc

Keep Your Lab Flowing:

Routine Maintenance Supplies for

Waters Alliance LCs

Replace Every Month

Bluestem Glass Solvent Filter

- Restek Bluestem glass solvent filter provides clean mobile phase to extend the life of columns and pump seals.
- 15 µm borosilicate glass frit sits lower than conventional glass filters to draw more mobile phase from each bottle.
- Blue filter stem allows instant visual confirmation of upright filter orientation.
- Connects to standard 1/8" OD (3.2 mm) PTFE tubing using your existing frit adaptor (also sold separately as cat.# 26392).

Prevent the particulates and microbial growth in your LC solvents from entering your instrument with the new Restek Bluestem glass solvent filter.

Glass Solvent Filter, 15 µm frit	5041-2168	ea.	26431
Frit Adaptor, PTFE	5062-8517	4-pk.	26392



Replace Every 3 Months

Vespel Rotor Seal

for Waters HPLC Systems: Rheodyne Style

Vespel Rotor Seal	1090, 7000, 7010, 7040, 7067	7010-039	ea.	25279
Vespel Rotor Seal	7125, 7126, 7725, 7725i, 9725	7125-047	ea.	25280

Seal Pack Rebuild Kit (without Seal Wash Tube Assembly)

for Waters HPLC Systems

Seal Pack Rebuild Kit, without Seal Wash Tube Assembly

Includes: wash tube seal; seal wash tube; PTFE washer; filter retainer; lower wash seal frit; needle wash frit; TFE washer; needle; injector seal (2); compression screw; ferrule	2690, 2695	WAT271019	kit	25495
---	------------	-----------	-----	-------

Head Plunger Seal Kits

for Waters HPLC Systems

Head Plunger Seals, Gold	2690, 2695, 2790, 2795	WAT270789, WAT270938	2-pk.	25387
Head Plunger Seals, Black	2690, 2695, 2790, 2795	WAT271066	2-pk.	25388
Head Plunger Seal Kit	ACQUITY, nano-ACQUITY	700002599	kit	26428



Replace Every 6 Months

Preventive Maintenance (PM) Kits for Waters HPLC & ACQUITY UPLC Systems

ACQUITY BSM Pump Kit

Includes: in-line filter assembly, stainless steel frit; tube assembly, transducer to check valve; tube assembly, SSV to in-line filter; primary check valve (2); wash seal, float flanged, (2); head plunger seal kit (2); sapphire plungers, (2); air filter, pump; air filter, pump handle; mixer assembly, 50 µL; accumulator check valve, double ball & seat (2); solvent bottle filter, stainless steel (7); pump O-ring, (PTFE)

ACQUITY UPLC
BSM


201000173

kit

25797



25797



25798

ACQUITY I2V BSM Pump Kit

Includes: tube assembly, SSV to I2v; transducer to check valve, tube assembly, I2v; wash seal, float flanged (2); head plunger seal kit (2); sapphire plungers (2); air filter, pump; air filter, pump handle; mixer assembly, 50 µL; filter frit cartridge, stainless steel; accumulator check valve, double ball & seat (2); check valve cartridge (2); solvent bottle filter, stainless steel (7); pump O-ring, (PTFE)

ACQUITY UPLC
I2V BSM

201000197

kit

25798



25143

ACQUITY H-Class QSM Pump Kit

Includes: solvent bottle filter, stainless steel (5); 20 micron frit holder assembly; tube assembly, transducer to check valve; wash seal, float flanged (2); head plunger seal kit (2); sapphire plungers (2); check valve, double ball & seat (2); I2 check valve cartridge; air filter, door; pump O-ring, (PTFE); mixer assembly, 100 µL

ACQUITY H-CLASS
QSM

201000233

kit

25799

2690/2695 Pump & Autosampler PM Kit

Includes: sapphire plungers (2); seal wash plunger seals (4); head plunger seals (4); wash tube seals (4); battery, sparge diffuser; filter insert; face seals (4); solvent reservoir 10 µm filters (4); 250 µL syringe, check valve cartridges (2); wash tube seal; seal wash tube; PTFE washer; filter retainer; lower wash seal frit; needle wash frit; TFE washer; needle assembly; gold injector seals (2); stainless steel ferrule; compression screw

Alliance 2690,
2695

WAT270944

kit

25143



25145

717 Autosampler PM Kit

Includes: seal pack assembly, tube assembly (0.020" ID), needle, needle compression screw, 0.062 stainless steel ferrule, pre-column filter assembly, filter insert, 250 µL WISP syringe

717 Autosampler

WAT052669

kit

25145

600 Pump PM Kit

Includes: PerformancePLUS cartridges (4); sparge diffusers (4); Super Seals (2); solvent reservoir 10 µm filters (4); sapphire plungers (2); reference valve button; valve disk spacer; valve disk; TFE ball plug; TFE seat; ruby ball; inlet tube body assembly manifold insert; insert seal; belleville washers (2); flat washer

600 Pump

WAT052675

kit

25144



25144



25146

616 Pump PM Kit

Includes: sapphire plungers (2); sparge diffusers (4); solvent reservoir 10 µm filters (4); plunger seals (2); plunger wash seals (2); O-rings (2); backup rings (2); check valve cartridges (4)

616 Pump

WAT052672

kit

25146

515 Pump PM Kit


Includes: PerformancePLUS check valves (4); sparge diffuser; solvent reservoir 10 µm filter; sapphire plungers (2); plunger seals (2); pivot inserts (2); pivot guides (2); washer (2); plunger springs (2); retaining rings (2)

515 Pump

WAT052587

kit

26519



26519



26430

1525 Pump PM Kit

Includes: sapphire plungers (4); check valve cartridges (8); plunger seals (4); solvent reservoir 10 µm filters (2); reference valve button; valve disk spacer; valve disk

1525 Pump

201000114

kit

26430

Check Valve Cartridges

for Waters HPLC Systems

Check Valve Cartridge	2690, 2695, 2790, 2795	WAT270941	2-pk.	25373
-----------------------	------------------------	-----------	-------	-------



25373

PerformancePLUS Check Valve Cartridge

for Waters HPLC Systems

PerformancePLUS Check Valve Cartridge	1515, 1525, 2690, 2695, 2795, 510, 515, 600, 610	700000254	2-pk.	25370
---------------------------------------	--	-----------	-------	-------



25370

Replace Every 12 Months

Sapphire Plungers

for Waters HPLC Systems

Sapphire Plunger	M6KA, 510, 590, 600, 610, LC Module 1	WAT025656, WAT069511	ea.	25381
Sapphire Plunger (Extended Flow)	1525EF, 510, 515, 590, 600, 610, LC Module 1	WAT060304	ea.	25382
Sapphire Plunger	616, 625, 626, 1525 Micro	WAT031788	ea.	25420
Sapphire Plunger	2690, 2695, 2790, 2795, 2796	WAT270488, WAT271067, WAT270959	ea.	25385



25382

Battery

for Waters HPLC Systems

Battery	2690, 2695, 717, 486, 484	WAT080443	ea.	25494
---------	---------------------------	-----------	-----	-------



25494

Don't Forget!

Lamps

for Waters Detectors

Lamp, UV, Long Life (2,000 hours)	ACQUITY PDA/TUV 2489/2998	201000281	ea.	25775
Deuterium Lamp, Long Life (2,000 hours)	486	WAT080678	ea.	25410
Deuterium Lamp, Standard (1,000 hours)	996, 2996	WAT052586	ea.	25408
Deuterium Lamp, Standard (1,000 hours)	2487, 2488	WAS081142	ea.	25409
Xenon Lamp (without holder or mirror)	470		ea.	25404
Xenon Lamp	474		ea.	25405



25409

This is just the beginning!

There are more RM supplies for your Waters, as well as parts for Agilent, PerkinElmer, Shimadzu, Thermo, and other LCs, at www.restek.com/LCacc